

HANDLING PIPETTE

MICROTOOL FOR TRANSFERRING EMBRYOS AND OOCYTES

The handling pipette is a microtool used for transferring oocytes and embryos from one dish to another.

FEATURES

- laboratory microtool for dish-to-dish transfer of oocytes and embryos
- will not scratch dishes
- available in three lengths: 45, 55 and 94 mm
- available in two packaging variations: single packed and packed by ten pieces



ORIGINAL PICTURE

GYNO INFO HANDLING PIPETTE

Length 45, 55 or 94 mm (\pm 3 mm)

Tube OD 1.25 mm (\pm 0.05 mm)

Tube ID 1.00 mm (\pm 0.05 mm)

Sterile sterilized by radiation

Packaging 250/shipping box
1250/shipping box (10/pouch)

Quality control MEA and LAL

Order number **#3145** ☒ (45 mm)
#3100 ☒ (55 mm)
#3125 ☒ (94 mm)
#3145 TP ☒ (45 mm, 10/pouch)
#3100 TP ☒ (55 mm, 10/pouch)
#3125 TP ☒ (94 mm, 10/pouch)



GMP Belgium Handling Pipette 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine